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InKina®

ENCRE DE CHINE 2 mL

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QUALITE

Téléphone : +39 328 033 5893

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AlaVia MEDS
16 avenue Rosa Parks
69009, Lyon – France
www.alaviameds.com

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EN

INSTRUCTIONS FOR USE

InKina[®]
INDIA INK 2 mL
 STERILE SURGERY MARKER FOR
 ENDOSCOPIC TATTOOING
 Sterile field ready

**NOT FOR INTRAVENOUS
 INJECTION OR ORAL ADMINISTRATION**

CE
0477

NAME OF MEDICAL DEVICE

InKina[®]

REFERENCE

MINS0002

TYPE OF MEDICAL DEVICE

Medical device of Class IIa.

QUALITATIVE AND QUANTITATIVE COMPOSITION

Black India ink.

USE

InKina[®] is a contrast dye indicated for the identification of pathological lesions in endoscopies of the gastrointestinal tract. The marker is used for the endoscopic tattooing of gastrointestinal mucosa lesions (colon lesions included), for subsequent surgical excision within the next 30 days.

SIDE EFFECTS AND CONTRAINDICATIONS

Clinical studies have shown that patients have not encountered any side effects or health disorders during the use of India ink. The use of the marker for tattooing the intraperitoneal area needs to be monitored to keep under control any undesirable side effects.

India ink is contraindicated in patients with known hypersensitivity towards the components of the marker.

WARNING : InKina[®] is not indicated for use in pregnant or breastfeeding patients or in pediatric patients.

INCOMPATIBILITY

Do not mix InKina[®] with medicines.

WARNINGS AND PRECAUTIONS FOR USE

**DO NOT INJECT InKina[®] INTRAVENOUSLY OR
 DO NOT ADMINISTER ORALLY**

INSTRUCTIONS FOR USE

InKina[®] is inoculated into the submucosa of the area affected by pathological lesions through a Luer-Lock needle during endoscopy.

InKina[®] can be diluted 1:10 in sterile saline solution.

The marked area is removed surgically.

- Before use, shake the marker vigorously in order to mix possible sediments formed.
- After stirring the suspension, connect the needle to the syringe. It is recommended to use a 25 Gauge sclerotherapy Luer-Lock needle with a length of 4 mm or less.
- Before injection, fill the catheter with the suspension.
- With the retracted needle, insert the catheter into the endoscope biopsy channel.
- When the needle is positioned appropriately, insert the tip with an inclination of 30°-45° into the submucosal layer of the colon.

WARNING : Do not insert the needle perpendicular to the colon, as it may puncture the serous membrane and insert InKina[®] directly into the peritoneal cavity.

- Inject about 0,5 – 0,75 mL of suspension at each injection site; if it is correctly performed, a black stain of 1-2 cm diameter is obtained.
- Observe if there is a leakage of suspension in the lumen of the colon. If this occurs, stop the injection and insert the needle deeper under the mucosa, paying close attention.
- A minimum of four injections is recommended, to be carried out circularly to the lesion to ensure that at least one of the markings can be observed by the surgeon from the colon serous membrane and, if another endoscopy is needed, the lesion can also be observed from within the lumen.
- It is recommended to use a maximum of 8 mL of suspension for each patient.

EXPIRY DATE

24 months referred to the marker in undamaged primary container.

Expiry date: Check the expiry date printed on the container.

The expiry date refers to the marker properly stored in an unopened container. Do not use the marker after the expiry date.

SPECIAL PRECAUTIONS FOR STORAGE

• Store in original container to protect the syringes from light. There is no special temperature storage conditions. Do not freeze.

- Do not use InKina[®] if the container is damaged.
- InKina[®] must be used for a single, uninterrupted administration and any residue must be discarded to prevent the risk of contamination due to loss of sterility.

KEEP OUT OF REACH AND SIGHT OF CHILDREN

CONTAINER

Primary : 1 prefilled syringe of 2 mL with a protective wrapping.

Secondary : one box of 5 syringes.

The graduated scale reported on the syringes are indicative and do not have a measurement value.

SPECIAL PRECAUTIONS FOR DISPOSAL AND OTHER HANDLING OF THE MARKER

The unused marker and the resulting waste must be disposed of in compliance with current local legislation.

SERIOUS INCIDENT REPORTING

The user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER

S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO

Via Marconi, 2 - 24069

Cenate Sotto (BG) – Italy

www.salfspa.it

VIGILANCE

Phone : +39 335 131 8126

E-mail : pharmacovigilance@salfspa.it

QUALITY

Phone : +39 328 033 5893

DISTRIBUTOR

AlaVia MEDS

16 avenue Rosa Parks

69009, Lyon – France

www.alaviameds.com

SYMBOLS

The symbols stated on the primary and secondary containers are in compliance with UNI CEI EN ISO 15223-1 rule (Symbols to be used with medical device labels, labelling and information to be supplied).

Symbols with relative descriptions are listed below:

	Medical Device Code
	Lot Number
	Expiry Date
	CE marking – Notified Body 0477
	Protect from direct sun light
	The marker does not contain latex
	Do not reuse / Single use only.
	Do not re-sterilize
	Pay attention to the instructions for use
	Sterilized using steam heat
	Do not use if the container is damaged
	Manufacturer
	Single sterile barrier system with protective packaging outside
	Caution
	Medical Device
	Throw into the appropriate container

REVISION DATE : February 2023

CXXXXX

Summary table of S.A.I.F. – AlaVia MEDS markers

NAME	FOR	USE	INSTRUCTIONS FOR USE	
 MethyBlue®	 BLUE BREAST SURGERY	In the visualization of sentinel lymph nodes to enable the excision of potentially malignant lymph nodes during surgery on the breast	Periareolar injection	2 to 5 mL of MethyBlue® undiluted (1%)
	UROLOGIC SURGERY	In the identification of the renal cavity during surgery or percutaneous nephrolithotomy	Intravesical injection	100 to 250 mL of MethyBlue® diluted (0,01%)
		Checking the tightness of urinary tract	Intravesical injection	25 to 250 mL of MethyBlue® diluted (0,01%)
	GASTRO-DIGESTIVE SURGERY	In the identification of the urinary or colonic fistula before surgical excision	Intravesical injection	10 to 100 mL of MethyBlue® diluted (0,01%)
		In the marking of the pilonidal sinus before surgical removal	Administer with a syringe directly into the pilonidal sinus	1 to 5 mL of MethyBlue® diluted (0,01%)
		In the preoperative marking of stoma prior to excision	Subcutaneous injection	2 to 5 mL of MethyBlue® diluted (0,01%)
Checking the tightness of colorectal sutures		Rectal administration	25 to 250 mL of MethyBlue® diluted (0,01%)	
 PURPLE LazuLymph®	GASTRO-DIGESTIVE SURGERY	Intraoperative sentinel lymph nodes mapping in colon cancer in early stage	Peritumoral injection	A total amount of less than 2 mL of LazuLymph® undiluted (2,5%)
	SKIN SURGERY	Preoperative identification of sentinel lymph nodes in melanoma	Peri-lesional intradermal injections 20 minutes prior to skin incision	0,5 to 1 mL of LazuLymph® undiluted (2,5%)
	BREAST SURGERY	Intraoperative sentinel lymph nodes mapping in patients with breast cancer	Subdermal, intradermal, periareolar, peritumoral or intratumoral injection	Subdermal : 0,2 to 0,4 mL of LazuLymph® undiluted (2,5%) Peritumoral, LazuLymph® undiluted (2,5%) : - 1 to 2 mL for tumors in the supero-external quadrant - 2 to 3 mL for tumors at the confluence of external or internal quadrants - 3 to 4 mL for tumors in lower quadrants.
 WHITE Inkina®	GASTRO-DIGESTIVE SURGERY	Identification of pathological lesions in endoscopies of the gastrointestinal tract.	Submucosa injection	Inkina® undiluted or diluted (1:10) 0,5 to 0,75 mL of suspension at each injection site Maximum of 8 mL of suspension for each patient

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Tableau récapitulatif des marqueurs de S.A.I.F. – AlaVia MEDS