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CXXXXX

INSTRUCTIONS FOR USE

LazuLymph®
PATENT BLUE V 2,5%
2 mL - 25mg/mL
MARKER FOR SURGERY
STERILE PYROGEN-FREE SOLUTION
 Sterile field ready

**NOT FOR INTRAVENOUS
 INJECTION**

CE
0477

NAME OF MEDICAL DEVICE

LazuLymph®

REFERENCE

MBPS2502

TYPE OF MEDICAL DEVICE

Sterile pyrogen-free solution – Class IIa.

QUALITATIVE AND QUANTITATIVE COMPOSITION

Patent blue V sodium salt 25 g – Sodium hydroxide q.s. pH 9,5 - Water for injection q.s. to 1000 mL.

USE

LazuLymph® is a vital dye suitable for the marking of sentinel lymph nodes in the breast cancer and in the intraoperative sentinel lymph nodes mapping in patients with colon cancer. It is also used in the preoperative lymphoscintigraphy for the biopsy of sentinel lymph nodes in the melanoma.

CONTRAINDICATIONS

Hypersensitivity to patent blue V 2,5% sodium salt.

WARNING: LazuLymph® should not be administered during pregnancy. Avoid the use during the lactation and in pediatric patients, because the related studies are not available.

SIDE EFFECTS

Patent blue V 2,5% sodium salt can cause what follows:

- hives, itch, rash;
- transient hypotension, bronchospasm / laryngospasm;
- rare cases of anaphylaxis, anaphylactic shock and angioneurotic edema.

WARNINGS AND PRECAUTIONS FOR USE

**LazuLymph® MUST NOT BE INJECTED
 INTRAVENOUSLY**

LazuLymph® can also cause severe allergic reactions. Before use, check the patient's allergies and intolerances.

Anesthesiologists and surgeons should carefully consider the onset of the first symptoms of an allergic reaction. They should also ensure the availability of personnel and equipment necessary to treat a possible allergic reaction, **supervising for at least 1 hour from the time of administration**, as there could be a possible delay of the onset of an adverse reaction. A blueish coloration of the surrounding skin may appear after the injection. It normally disappears within 24 - 48 hours. In case of an intradermal injection the coloration of the skin may persist for a longer period of time.

The dye is eliminated in 24 - 48 hours after the injection via the urine which are consequently strongly colored, and also via the bile.

INSTRUCTIONS FOR USE

To be used with a Luer Lock connection needle:

Intraoperative sentinel lymph nodes mapping in patients with breast cancer

The vital dye may be injected by subdermal, intradermal, periareolar, peritumoral or intratumoral route, although this last route is rarely used, because it requires the administration of high volumes of dye, given the paucity of lymphatics within the tumor mass.

The amount to be used depends on the route of administration:

- By subdermal route you should inject (approximately) 0,2 – 0,4 mL.
- By peritumoral route you should inject larger volumes in two or more peritumoral points: (approximately) 1 - 2 mL for tumors in the supero-external quadrant; (approximately) 2 - 3 mL for tumors at the confluence of external or internal quadrants; (approximately) 3 - 4 mL for tumors in lower quadrants.

Intraoperative sentinel lymph nodes mapping in colon cancer in early stage

A total amount of less than 2 mL of vital dye LazuLymph® is injected into the peritumoral region, at the four cardinal points, in the subserosal area to allow the identification of the first 4 sentinel nodes.

Preoperative identification of sentinel lymph nodes in melanoma

Peri-lesional intradermal injections of (approximately) 0,5 - 1 mL of LazuLymph® are made 20 minutes prior to skin incision.

EXPIRY DATE

24 months referred to the marker in undamaged container.

Expiry date: Check the expiry date printed on the container. The expiry date refers to the marker properly stored in an unopened primary container.

Do not use the marker after the expiry date.

SPECIAL PRECAUTIONS FOR STORAGE

- Store in original container to protect the syringes from light. There are no special temperature storage conditions. Do not freeze.
- Do not use LazuLymph® if the container is damaged.
- LazuLymph® must be used for a single, uninterrupted administration and any residue must be discarded to prevent the risk of contamination due to loss of sterility.

KEEP OUT OF REACH AND SIGHT OF CHILDREN**CONTAINER**

Primary: 1 prefilled syringe of 2 mL with a protective wrapping.

Secondary: one box of 5 syringes.

The graduated scale reported on the syringes are indicative and do not have a measurement value.

SPECIAL PRECAUTIONS FOR DISPOSAL AND OTHER HANDLING OF THE MARKER

The unused marker and the resulting waste must be disposed of in compliance with current local legislation.

SERIOUS INCIDENT REPORTING

The user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER

S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO

Via Marconi, 2 - 24069

Cenate Sotto (BG) – Italy

www.salfspa.it

**VIGILANCE**

Phone: +39 335 131 8126

E-mail: pharmacovigilance@salfspa.it

QUALITY

Phone: +39 328 033 5893

DISTRIBUTOR

AlaVia MEDS

16 avenue Rosa Parks

69009, Lyon – France

www.alaviameds.com

**SYMBOLS**

The symbols stated on the primary and secondary containers are in compliance with UNI CEI EN ISO 15223-1 rule (Symbols to be used with medical device labels, labelling and information to be supplied). Symbols with relative descriptions are listed below:

	Medical Device Code
	Lot Number
	Expiry Date
	CE marking – Notified Body 0477
	Protect from direct sun light
	The marker does not contain latex
	Do not reuse / Single use only.
	Do not re-sterilize
	Pay attention to the instructions for use before using the marker
	Sterilized using steam heat
	Do not use if the container is damaged
	Manufacturer
	Single sterile barrier system with protective packaging outside
	Caution
	Medical Device
	Throw into the appropriate container

REVISION DATE: February 2023

CXXXXX

Summary table of S.A.L.F. – AlaVia MEDS markers

NAME	FOR	USE	INSTRUCTIONS FOR USE	
 MethyBlue®	 BREAST SURGERY	In the visualization of sentinel lymph nodes to enable the excision of potentially malignant lymph nodes during surgery on the breast	Periareolar injection	2 to 5 mL of MethyBlue® undiluted (1%)
		UROLOGIC SURGERY	In the identification of the renal cavity during surgery or percutaneous nephrolithotomy	Intravesical injection
	Checking the tightness of urinary tract		Intravesical injection	25 to 250 mL of MethyBlue® diluted (0,01%)
	GASTRO-DIGESTIVE SURGERY	In the identification of the urinary or colonic fistula before surgical excision	Intravesical injection	10 to 100 mL of MethyBlue® diluted (0,01%)
		In the marking of the pilonidal sinus before surgical removal	Administer with a syringe directly into the pilonidal sinus	1 to 5 mL of MethyBlue® diluted (0,01%)
		In the preoperative marking of stoma prior to excision	Subcutaneous injection	2 to 5 mL of MethyBlue® diluted (0,01%)
Checking the tightness of colorectal sutures	Rectal administration	25 to 250 mL of MethyBlue® diluted (0,01%)		
 Lazulymph®	GASTRO-DIGESTIVE SURGERY	Intraoperative sentinel lymph nodes mapping in colon cancer in early stage	Peritumoral injection	A total amount of less than 2mL of Lazulymph® undiluted (2,5%)
	SKIN SURGERY	Preoperative identification of sentinel lymph nodes in melanoma	Peri-lesional intradermal injections 20 minutes prior to skin incision	0,5 to 1 mL of Lazulymph® undiluted (2,5%)
	BREAST SURGERY	Intraoperative sentinel lymph nodes mapping in patients with breast cancer	Subdermal, intradermal, periareolar, peritumoral or intratumoral injection	Subdermal : 0,2 to 0,4 mL of Lazulymph® undiluted (2,5%) Peritumoral, Lazulymph® undiluted (2,5%) : - 1 to 2 mL for tumors in the supero-external quadrant - 2 to 3 mL for tumors at the confluence of external or internal quadrants - 3 to 4 mL for tumors in lower quadrants.
 Inkina®	GASTRO-DIGESTIVE SURGERY	Identification of pathological lesions in endoscopies of the gastrointestinal tract.	Submucosa injection	Inkina® undiluted or diluted (1:10) 0,5 to 0,75 mL of suspension at each injection site Maximum of 8 mL of suspension for each patient

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Tableau récapitulatif des marqueurs de S.A.L.F. – AlaVia MEDS