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QUALITE

Téléphone : +39 328 033 5893

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AlaVia MEDS
16 avenue Rosa Parks
69009, Lyon – France
www.alaviameds.com

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INSTRUCTIONS FOR USE
MethyBlue®
METHYLENE BLUE 1%
5 mL - 10mg/mL
STERILE SURGICAL DYE
Sterile field ready

**NOT FOR INTRAVENOUS,
 INTRATHECAL OR
 INTRA-AMNIOTIC INJECTION**

CE
0477

NAME OF MEDICAL DEVICE
 MethyBlue®

REFERENCE
 MBMS0105

TYPE OF MEDICAL DEVICE
 Sterile pyrogen-free solution – Class IIa.

QUALITATIVE AND QUANTITATIVE COMPOSITION

Methylthionium chloride trihydrate (methylene blue trihydrate) 10 g - Hydrochloric acid or sodium hydroxide to adjust pH as needed - Water for injection q.s. to 1000 mL.

USE

MethyBlue® was designed to be used in surgery as a way of marking tissues and operative findings and to check the tightness of urinary tract and colorectal sutures and to evaluate tubal patency.

CONTRAINDICATIONS

Known hypersensitivity to methylthionium chloride trihydrate (methylene blue).

WARNING: Do not use MethyBlue® during pregnancy.

WARNINGS AND PRECAUTIONS FOR USE

MethyBlue® IS NOT FOR INTRAVENOUS, INTRATHECAL OR INTRA-AMNIOTIC INJECTION.

For the described indications, MethyBlue® is not a specific marker. It is therefore **not to be used for diagnostic purposes** in humans, for example in screening for precancerous lesions.

Adverse effects vary according to the use of the medical device.

When MethyBlue® is administered by subcutaneous injection, it can cause vasoconstriction or skin necrosis. In case of contact with mucous membranes in the eyes, flush thoroughly with plenty of water. The use of MethyBlue® may cause temporary colouring of urine and faeces and integuments. The following adverse effects were also reported:

- hypersensitivity and anaphylactic reaction, pulmonary edema;
- pre- or postoperative methemoglobinemia.

INSTRUCTIONS FOR USE

According to its intended use, the device MethyBlue can be connected with either a needle or a catheter. Both of them must be equipped with a standard Luer Lock connection to ensure compatibility.

MARKING OF TISSUES AND OPERATIVE FINDINGS

In the visualization of sentinel lymph nodes to enable the excision of potentially malignant lymph nodes during breast surgery

Administer by periareolar injection a quantity between (approximately) 2 mL and 5 mL MethyBlue® 1%.

For all other indications listed below, use a methylene blue solution with a concentration of 0,01%. Therefore, before use, dilute 1 part of MethyBlue® with 100 parts of sterile sodium chloride 0,9%.

In the identification of the renal cavity during surgery or percutaneous nephrolithotomy Administer intravesically a quantity of diluted MethyBlue® between 100 mL and 250 mL.

In the identification of the urinary or colonic fistula before surgical excision Administer intravesically a quantity of diluted MethyBlue®, generally of between 10 mL and 100 mL.

In the marking of the pilonidal sinus before surgical removal Administer a quantity of diluted MethyBlue®, generally between 1 mL and 5 mL according to the size of the cyst directly into the pilonidal sinus with a syringe.

In the preoperative marking of stoma prior to excision

The volume of diluted MethyBlue® to inject subcutaneously is at the physician's discretion (generally between 2 mL and 5 mL).

CHECKING THE TIGHTNESS OF URINARY TRACT AND COLORECTAL SUTURES

Use a methylene blue solution with a concentration of 0,01%. Therefore, before use, dilute 1 part of MethyBlue® with 100 parts of sterile sodium chloride 0,9%.

Administer intravesically or rectally a quantity of diluted MethyBlue® between 25 mL and 250 mL according to the position of the suture. If no blue colouring is detected in the suture area, then the suture is proven to be properly sealed.

EVALUATION OF TUBAL PATENCY:

Administer the diluted solution (approximately 10 mL) via the intrauterine route. The flow of diluted methylene blue through the Fallopian tubes confirms tubal patency.

EXPIRY DATE

36 months in an unopened primary container. Expiry date: Check the expiry date printed on the container.

The expiry date refers to the marker properly stored in an unopened primary container. Do not use the marker after the expiry date.

SPECIAL PRECAUTIONS FOR STORAGE

- Store in original container to protect the syringes from light. There is no special temperature storage conditions. Do not freeze.
- Do not use MethyBlue® if the primary container is damaged.
- MethyBlue® must be used for a single, uninterrupted administration and any residue must be discarded to prevent the risk of contamination due to loss of sterility.

KEEP OUT OF REACH AND SIGHT OF CHILDREN

CONTAINER

Primary: 1 prefilled syringe of 5 mL with a protective wrapping.
 Secondary: one box of 5 syringes.
 The graduated scale reported on the syringes are indicative and do not have a measurement value.

SPECIAL PRECAUTIONS FOR DISPOSAL AND OTHER HANDLING OF THE MARKER

The unused marker and the resulting waste must be disposed of in compliance with current local legislation.

SERIOUS INCIDENT REPORTING

The user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER

S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO
 Via Marconi, 2 - 24069
 Cenate Sotto (BG) – Italy
www.salfspa.it

VIGILANCE

Phone : +39 335 131 8126
 E-mail : pharmacovigilance@salfspa.it

QUALITY

Phone: +39 328 033 5893

DISTRIBUTOR

AlaVia MEDS
 16 avenue Rosa Parks
 69009, Lyon – France
www.alaviameds.com

SYMBOLS

The symbols stated on the primary and secondary containers are in compliance with UNI CEI EN ISO 15223-1 rule (Symbols to be used with medical device labels, labelling and information to be supplied).

Symbols with relative description are listed below:

	Medical Device Code
	Lot Number
	Expiry Date
	CE marking – Notified Body 0477
	Protect from direct sun light
	The marker does not contain latex
	Do not reuse / Single use only.
	Do not re-sterilize
	Pay attention to the instructions for use before using the marker
	Sterilized using steam heat
	Do not use if the container is damaged
	Manufacturer
	Single sterile barrier system with protective packaging outside
	Caution
	Medical Device
	Throw into the appropriate container

REVISION DATE : June 2024

Summary table of S.A.L.F. – AlaVia MEDS markers				
NAME	FOR	USE	INSTRUCTIONS FOR USE	
	BREAST SURGERY	In the visualization of sentinel lymph nodes to enable the excision of potentially malignant lymph nodes during surgery on the breast	Periareolar injection	2 to 5 mL of MethyBlue® undiluted (1%)
		In the evaluation of tubal patency	Intrauterine route	Approximately 10 mL of MethyBlue® diluted (0,01%)
	UROLOGIC SURGERY	In the identification of the renal cavity during surgery or percutaneous nephrolithotomy	Intravesical injection	100 to 250 mL of MethyBlue® diluted (0,01%)
		Checking the tightness of urinary tract	Intravesical injection	25 to 250 mL of MethyBlue® diluted (0,01%)
	GASTRO-DIGESTIVE SURGERY	In the identification of the urinary or colonic fistula before surgical excision	Intravesical injection	10 to 100 mL of MethyBlue® diluted (0,01%)
		In the marking of the pilonidal sinus before surgical removal	Administer with a syringe directly into the pilonidal sinus	1 to 5 mL of MethyBlue® diluted (0,01%)
In the preoperative marking of stoma prior to excision		Subcutaneous injection	2 to 5 mL of MethyBlue® diluted (0,01%)	
	Checking the tightness of colorectal sutures	Rectal administration	25 to 250 mL of MethyBlue® diluted (0,01%)	
	GASTRO-DIGESTIVE SURGERY	Intraoperative sentinel lymph nodes mapping in colon cancer in early stage	Peritumoral injection	A total amount of less than 2 mL of Lazulymph® undiluted (2,5%)
	SKIN SURGERY	Preoperative identification of sentinel lymph nodes in melanoma	Peri-lesional intradermal injections 20 minutes prior to skin incision	0,5 to 1 mL of Lazulymph® undiluted (2,5%)
	BREAST SURGERY	Intraoperative sentinel lymph nodes mapping in patients with breast cancer	Subdermal, intradermal, periareolar, peritumoral or intratumoral injection	Subdermal : 0,2 to 0,4 mL of Lazulymph® undiluted (2,5%) Peritumoral, Lazulymph® undiluted (2,5%) : - 1 to 2 mL for tumors in the supero-external quadrant - 2 to 3 mL for tumors at the confluence of external or internal quadrants - 3 to 4 mL for tumors in lower quadrants.
	GASTRO-DIGESTIVE SURGERY	Identification of pathological lesions in endoscopies of the gastrointestinal tract.	Submucosa injection	Inkina® undiluted or diluted (1:10) 0,5 to 0,75 mL of suspension at each injection site Maximum of 8 mL of suspension for each patient

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Tableau récapitulatif des marqueurs de S.A.L.F. – AlaVia MEDS